


**Compact  
Powerful  
Precise**



# ColorFlex<sup>®</sup> L2

*Level UP your Color Confidence with the new  
ColorFlex L2 Spectrophotometer.*

**As simple as you want, as powerful as you need.**

Elevate your color accuracy and user experience like never before. Designed to streamline and expedite your quality control processes, ColorFlex L2 Spectrophotometer gives you unsurpassed confidence and trust in your outcomes. Compact and ergonomic build requires minimal space while offering enhanced durability and reliability.


# Explore the Next-Gen ColorFlex®...L2



Precise color measurement for everything from opaque solids, liquids, powders, granules and pellets to translucent solids and liquids. We have upgraded the most relied upon 45°/0° spectrophotometer in the world, ColorFlex EZ, to the new ColorFlex L2 with the most advanced color technology.



**Explore all the additional features and benefits below.**



## 45°a:0° Annular Illuminated Optical Geometry

**Superior Consistency and Color Accuracy:** Enhancing the efficiency of Quality Control workflows and delivery of superior quality products repeatedly.



## Compact with Sealed Spillproof Case

**Space-Saving and Easy Maintenance:** Occupies minimal space and is designed to withstand the demands of busy lab environments.

## Powerful Standalone Functionality

**No computer required:** With huge onboard processing and storage capacity, ColorFlex L2 can function as an independent color workstation. Connect monitor, keyboard, and mouse if desired.




Sliced Cheese      SPICES      flour      POWDERED SPORTS DRINKS      Yogurt  
Calcium Carbonate      Puree      mustards      barium sulfate  
BABY FOOD      Coffee - Green beans, Roasted & Ground      CEMENT  
Baking Powders      Citrus Juice      Herbs      Cocoa Powder



## Modern Tablet User Interface with Touchscreen & Color Display

**Durable and Efficient:** Facilitates accurate color measurement in seconds, with a robust industrial glass design suitable for harsh environments.




## New & Improved Color Software

**Essentials 2.0:** Take precise and reliable color measurement in seconds with next-gen user-friendliness and seamlessly analyze multiple tristimulus and spectral data sets. It's that powerful.


## Self-Guided and Intuitive Out-of-Box Startup

**Time-Saving and User-Friendly:** Ready to use straight out of the box with an integrated wizard for quick training, reducing operator time and easing concerns related to staff training and turnover.



## Built-in Camera for Sample Viewing and Screen Recording

**Precision and Communication:** Ensures accurate sample positioning to eliminate measurement discrepancies and allows easy sharing of recordings for swift quality control decisions.



## Built-in Advanced Diagnostics

**On-Board Diagnostic Health Check:** Ensures the HunterLab ColorFlex L2 is always operating at peak performance, continuing to build on the legendary history of reliability of the ColorFlex family of instruments.

Powdered Sugar	BAKING CHIPS	Concrete	Powdered Functional Beverages	Fertilizers
Magnesium Sulfate	tea drinks (powder)	Starch	shortening	WHITE SUGAR
RAW GRAINS	Brown Sugar	powdered cheese	SEASONINGS	Orange Juice



## Advanced Data Export and Communication

**Full suite of data communication options:** Features USB ports, HDMI, and support for mouse and keyboard, enabling easy data export/printing.



## Connectivity with LIMS and SPC systems

**Simplifies QC workflows:** Aids in trend analysis and workflow simplification across manufacturing sites.



# CONTACT US

We've worked with thousands of companies, but you'll feel like our one and only. If color accuracy is crucial to your business, you've come to the right place—with the right people. Let's get in touch. Ask a question, Request a quote or just say hello!

## OUR OFFICES

### HunterLab

11491 Sunset Hills Road  
Reston, VA 20190-5280  
Website: [www.hunterlab.com](http://www.hunterlab.com)  
Email: [sales@hunterlab.com](mailto:sales@hunterlab.com)  
Ph: 703-471-6870  
Fax: 703-471-4237

### Corporate Headquarters

### HunterLab Europe

Christian Jansen, Managing Director  
Dr.-August-Einsele-Ring 15  
82418 Murnau, Germany  
Website: [www.hunterlab.de](http://www.hunterlab.de)  
Email: [Christian.Jansen@hunterlab.com](mailto:Christian.Jansen@hunterlab.com)  
Ph: +49 8841 9464

### European Office

### HunterLab Asia-Pacific

Leo Cheung, Regional Manager  
Flat E, 13F Leahander Centre  
28 Wang Wo Tsai Street Tsuen Wan  
New Territories Hong Kong  
Website: [www.hunterlabcolor.cn](http://www.hunterlabcolor.cn)  
Email: [sales-asiapacific@hunterlab.com](mailto:sales-asiapacific@hunterlab.com)  
Ph: +852-2406 1723

### Asia-Pacific Office

For more information, please contact HunterLab at 703-471-6870, [sales@hunterlab.com](mailto:sales@hunterlab.com) or visit [www.hunterlab.com](http://www.hunterlab.com)